



PLASTIC BOX

TRANSPORT BOX



RIGID BLISTER

SHELF BOX

TRANSPORT BOX

STERILE

4891

Currettes

Curretas

Currettes


Kùrettes

Currettes




Currettes

Curretas


Currettes



4891-001/501



4891-002/502



GREY CURETTES

- High performance polymer.
- Used for scraping, debriding and/or cleaning of skin and soft tissues and wounds.
- Two models:
 - Double Curette.
 - Curette with scraper.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton
 - Sterile: Individual peel-open pack + shelf box + transport carton.

CURETAS CINZA

- Polímero de alta resistência.
- Utilizado para raspagem, desbridamento e/ou limpeza de pele, tecidos moles e feridas.
- Dois modelos:
 - Cureta dupla.
 - Cureta com raspador.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizadas: Caixa de plástico + caixa de transporte.
 - Esterilizadas: Embaladas individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-001	Curette with scraper /Cureta c/raspador	14cm	800	1.600
4891-002	Double Curette/Cureta dupla	14cm	800	1.600

Sterile • Esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-501	Curette with scraper /Cureta c/raspador	14cm	1	30	120
4891-502	Double Curette/Cureta dupla	14cm	1	30	120