



Adson forceps

Pinça de Adson

Pincas Adson

Adson-Pinzette

Adson-pincet

Pinza Adson




Fórceps Adson

Adson Pincett



NON-TOOTHED
SEM DENTE
NON DENTÉE
NITCHZAHN
ZONDER TANDEN
NON DENTATE
SIN DIENTE
ICKE-TANDAD

1X2 TEETH
1X2 DENTE
1X2 GRIFFES
1X2 ZAHNE
1X2 TANDEN
1X2 DENTATO
1X2 DIENTES
KLO 1X2



ADSON TISSUE FORCEPS GREY

- High performance polymer.
- Intended for grasping tissues and/or holding other devices.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

PINÇA DE ADSON CINZA

- Polímero de alta performance.
- Indicadas para apreensão de tecidos e/ou manusear outros dispositivos.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizadas: Caixa de plástico + caixa de transporte.
 - Esterilizadas: Embalados individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-003	Non-toothed / Sem dentes	12cm	208	416
4891-004	1x2 teeth / 1x2 dentes	12cm	208	416

Sterile • Esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-504	Non-toothed / Sem dentes	12cm	1	15	60
4891-505	1x2 teeth / 1x2 dentes	12cm	1	15	60