

4891**ZILLION[®]**
BLACK

PLASTIC BOX

TRANSPORT BOX



RIGID BLISTER

SHELF BOX

TRANSPORT BOX

STERILE



Peaen Forceps

Pinça Pean

Pince Pean

Pean-Klemme

Peanklem

Pinza Pean

Fórceps de Pean

Pean peang



STRAIGHT
RECTA
DROIT
GERADE
RECHT
PUNTA
PUNTA
RAK



CURVED
CURVA
COURBE
GEBOGEN
GEBOGEN
CURVO
CURVA
BÓJD



PEAN FORCEPS GREY

- High performance polymer.
- Intended for hemostasis/clamp vessels and/or grasping tissues.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

PINÇA PEAN CINZA

- Polímero de alta performance.
- Indicadas para hemostase/clampagem de vasos e/ou prensão de tecidos.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizadas: Caixa de plástico + caixa de transporte.
 - Esterilizadas: Embaladas individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-009	Straight /Recta	14,5cm	144	288
4891-010	Curved/Curva	14,5cm	144	288

Sterile • Esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-510	Straight /Recta	14,5cm	1	33	132
4891-511	Curved/Curva	14,5cm	1	24	96

Bastos Viegas,s.a.

Cat, BV10, Rev 0 (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com