

ZILLION[®]
BLACK



PLASTIC BOX



TRANSPORT BOX



RIGID BLISTER



SHELF BOX



TRANSPORT BOX



STERILE

STRAIGHT
RECTA
DROIT
GERADE
RECHT
PUNTA
PUNTA
RAK

4891

Spencer-Wells forceps

Pinça Spencer-Wells

Pince Spencer-Wells

Spencer-Wells-Klemme

Spencer-wellsklem

Pinza Spencer-wells

Fórceps Spencer-wells

Spencer-wells peang



CURVED
CURVA
COURBE
GEBOGEN
GEBOGEN
CURVO
CURVA
BÓJD



SPENCER-WELLS FORCEPS GREY

- High performance polymer.
- Intended for hemostasis/clamp vessels and/or grasping tissues.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

PINÇA SPENCER-WELLS CINZA

- Polímero de alta performance.
- Indicadas para hemostase/clampagem de vasos e/ou prensão de tecidos.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizadas: Caixa de plástico + caixa de transporte.
 - Esterilizadas: Embaladas individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-015	Straight /Recta	18,5cm	100	200
4891-016	Curved/Curva	18,5cm	100	200

Sterile • Esterilizadas

REF	Description Descrição	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-516	Straight /Recta	18,5cm	1	33	132
4891-517	Curved/Curva	18,5cm	1	26	104

Bastos Viegas,s.a.

Cat. BV10, Rev O (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com