

4891

ZILLION[®]
BLACK

STERILE



STRAIGHT
RECTA
DROIT
GERADE
RECHT
PUNTA
PUNTA
RAK

Needle holder Mayo-Hegar

Porta-agulhas Mayo-Hegar

Porte aiguille Mayo-Hegar


Mayo hegar nadelhalter

Naaldvoerder Mayo Hegar

Porta aghi Mayo-Hegar

Soporte para agujas Mayo-Hegar

Nålförare Mayo Hegar



NEEDLE HOLDER MAYO-HEGAR GREY

- High performance polymer.
- Intended to hold and guide suture needles securely for suturing.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

PORTA-AGULHAS MAYO-HEGAR CINZA

- Polímero de alta performance.
- Destina-se a segurar e orientar as agulhas de forma segura durante a sutura.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizado: Caixa de plástico + caixa de transporte.
 - Esterilizado: Embalado individualmente em peel-open pack + caixa de cartolina + caixa de transporte

Non sterile • Não esterilizados

REF	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-019	15cm	120	240

Sterile • Esterilizados

REF	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-520	15cm	1	33	132

Bastos Viegas,s.a.

Cat, BV10, Rev O (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com