

4891

ZILLION[®]
BLACK

PLASTIC BOX

TRANSPORT BOX



RIGID BLISTER

SHELF BOX

TRANSPORT BOX



STERILE

Skin staple remover

Removedor de agrafos

Ôte-agrafes cutanées


Hautklammerentferner

Huidnietjesverwijderaar

Rimuovi graffe

Quitagrapas

Agraffborttagare



SKIN STAPLE REMOVER GREY

- High performance polymer.
- Intended to provide removal of surgical skin staples.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

REMOVEDOR DE AGRAFOS CINZA

- Polímero de alta performance.
- Destina-se a remover agrafos cirúrgicos da pele.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizado: Caixa de plástico + caixa de transporte.
 - Esterilizado: Embalados individualmente em peel-open pack + caixa de cartolina + caixa

Non sterile • Não esterilizado

REF	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-020	12,5cm	160	320

Sterile • Esterilizado

REF	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-521	12,5cm	1	28	112

Bastos Viegas,s.a.

Cat, BV10, Rev 0 (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com