

ZILLION[®]
BLACK

Sponge forceps

Pinças de gaze

Porte tampons

Tupferklemme

Depperklem

Pinza ad Anelli

Pinza para apósitos

Tvättpeang



PLASTIC BOX

TRANSPORT BOX



RIGID BLISTER

SHELF BOX

TRANSPORT BOX



STERILE



Sponge Forceps

- Produced with a High performance polymer.
- Intended to hold swabs for absorption during surgery or wound care procedures, and for skin disinfection procedures.
- Sterilize by moist heat at 121°C or 134°C*, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton
 - Sterile: Individual peel-open pack + shelf box + transport carton.

PINÇAS DE GAZE

- Fabricadas num polímero de alta performance.
- Indicadas para segurar compressas para absorção durante procedimentos cirúrgicos ou tratamento de feridas e para procedimentos de desinfecção da pele.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizados: Caixa de plástico + caixa de transporte.
 - Esterilizados: Embalados individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizadas

REF	Length Comprimento	Description Descrição	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-030	26,5cm	Curved/Curva	45	90
4891-033	20cm	Straight/Recta	80	160
4891-034	20cm	Curved/Curva	80	160

Sterile • Esterilizadas

REF	Length Comprimento	Description Descrição	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-531	26,5cm	Curved/Curva	1	9	81
4891-534	20cm	Straight/Recta	1	30	120
4891-535	20cm	Curved/Curva	1	7	28