

4891

ZILLION[®]
BLACK

PLASTIC BOX

TRANSPORT BOX



RIGID BLISTER

SHELF BOX

TRANSPORT BOX



STERILE

Roux retractor

Afastador Roux

Écarteur de Roux

Roux-Haken

Roux-retractor

Divaricatore Roux

Separador de Roux

Roux-sårhake



ROUX RETRACTOR GREY

- High performance polymer.
- Intended to hold back tissues and organs in order to expose the surgical field.
- Sterilizable by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton.
 - Sterile: Individual peel-open pack + shelf box + transport carton.

AFASTADOR ROUX CINZA

- Polímero de alta performance.
- Indicado para segurar tecidos e órgãos, a fim de expor o campo cirúrgico.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizados: Caixa de plástico + caixa de transporte.
 - Esterilizados: Embalados individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizado

REF	Length Comprimento	Pcs./Plastic box Uni./Cx. plástico	Pcs./Tr. carton Uni./Cx. transporte
4891-031	17cm	90	180

Sterile • Esterilizado

REF	Length Comprimento	Pcs./Pack Uni./Pacote	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr. carton Uni./Cx. transporte
4891-532	17cm	1	7	28

Bastos Viegas,s.a.

Cat, BV10, Rev 0 (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com