


4891

ZILLION[®]
BLACK

STERILE



STRAIGHT
RECTA
DROIT
GERADE
RECHT
PUNTA
PUNTA
RAK



Surgical scissor

Tesoura cirúrgica

Ciseaux chirurgicaux

Chirurgische Schere

Chirurgische schaar

Forbice chirurgica

Tijeras quirúrgicas

Operations sax

SURGICAL SCISSOR

- Produced with a High performance polymer.
- To cut devices during dressing or surgical procedures.
- Sterilize by moist heat at 121°C or 134°C, gas (Ethylene Oxide).
- Packing:
 - Non sterile: Plastic box + transport carton
 - Sterile: Individual peel-open pack + shelf box + transport carton.

TESOURA CIRÚRGICA

- Fabricadas num polímero de alta performance.
- Indicadas para cortar dispositivos em procedimentos de penso ou cirúrgicos.
- Esterilizar por calor húmido a 121°C ou a 134°C, gás (Óxido de Etileno).
- Empacotamento:
 - Não esterilizados: Caixa de plástico + caixa de transporte.
 - Esterilizados: Embalados individualmente em peel-open pack + caixa de cartolina + caixa de transporte.

Non sterile • Não esterilizados

REF	Description Descrição	Length Comprimento	Pcs./Plastic box Unl./Cx. plástico	Pcs./Tr. carton Unl./Cx. transporte
4891-200	Straight/Recta SH/BL	14,5cm	120	240
4891-201	Straight/Recta SH/SH	14,5cm	120	240
4891-202	Straight/Recta BL/BL	14,5cm	120	240

Sterile • Esterilizados

REF	Length Comprimento	Length Comprimento	Pcs./Pack Unl./Pacote	Pcs./Sh.box Unl./Cx.cartolina	Pcs./Tr. carton Unl./Cx. transporte
4891-700	Straight/Recta SH/BL	14,5cm	1	24	96
4891-701	Straight/Recta SH/SH	14,5cm	1	24	96
4891-702	Straight/Recta BL/BL	14,5cm	1	24	96

Bastos Viegas,s.a.

Cat, BV10, Rev O (10-2021)

Av. da Fábrica 298, 4560-164 Guilhufe, Penafiel, Portugal . Tel/Phone: (+351) 255 729 500 - geral@bastosviegas.com - bastosviegas.com